# **Electroencephalography - An Introduction** Klaus Lehnertz (Summer Term) mathing and a state of the second www.www.www

## Nomenclature

## **Electroencephalography** (EEG)

*"the recording of electrical activity along the scalp produced by the firing of neurons within the brain"* 

in clinical context:

*"the recording of the brain's spontaneous electrical activity over a short (sometimes long) period of time as recorded from multiple electrodes placed on the scalp"* 

## **Some Historical Milestones**

- 1791 Luigi Galvani: electrical stimulation of frog nerves
- 1849 Hermann von Helmholtz: speed of frog nerve impulses

#### 1850 - Emil Du Bois-Reymond: nerve galvanometer

- 1868 Julius Bernstein: time course of action potential
- 1870 Eduard Hitzig and Gustav Fritsch discover cortical motor area of dog using electrical stimulation
- 1874 Roberts Bartholow: electrical stimulation of human cortical tissue
- 1875 Richard Caton: recordings of electrical activity from an animal brain
- 1906 Sir Charles Scott Sherrington publishes

*The Integrative Action of the Nervous System* that describes the synapse and motor cortex


1913 - Edgar Douglas Adrian: all-or-none principle in nerve

#### 1929 - Hans Berger: first human electroencephalogram

- 1932 Jan Friedrich Toennies: multichannel ink-writing EEG machine
- 1932 Nobel Prize to E. D. Adrian and C. S. Sherrington for work on the function of neurons
- 1932 Jan Friedrich Toennies and Brian Matthews: differential amplifier
- 1991 Nobel Prize to Erwin Neher and Bert Sakmann for work on the function of single ion channels




## Some Useful Facts about the Brain




#### Cortex:

thickness: 1.5 - 4.5 mm total area: 2200 cm<sup>2</sup> (2/3 in sulci) # neurons: ~  $10^{10}$  (75 % excitatory) # synapses / neuron: ~  $10^3 - 10^4$ length of connections: ~  $10^7$ -  $10^9$  m (~ 2.5 x distance earth - moon) # ion channels / neuron: ~  $10^2 - 10^3$ exchange of information: electromagnetic, chemical neurotransmitter & other active substances: ~ 50 glia cells: x3 # neurons

## Cyto-architectural Map of Human Cerebral Cortex




#### **Brodman areas**




## Important organizational features of cortex


#### horizontal lamination




#### vertical columnation



## Inputs-outputs to Different Layers are Different



## EEG is Generated by Large Pyramidal-shaped Neurons in Layers II, III, IV




#### Important features:

- oriented vertical to cortical surface
- inhibitory and excitatory inputs are spatially segregated over the surface of these neurons

soma: inhibitory inputs only dendrites: excitatory and inhibitory inputs (excitatory : inhibitory ~ 6.5 : 1)

inhibitory inputs to dendrites and soma not generally the same

## How Charge Separation arises in Cortical Neurons?



#### Assumption:


excitation is increased in dendrites (other choices are possible)

Generators of EEG are excitatory and inhibitory post-synaptic potentials (not action potentials !!)


#### **Convention:**

current flow is the direction that the positive ions flow

## Current flow due to EPSP and IPSP



## Field potentials and EEG



neuron as *current dipole* (current flow between soma and apical dendrites)

de- or hyperpolarization of postsynaptic membrane

cell body differentially polarized

polarization travels to cell body

-> field potential

potential difference can be measured between electrodes located at different iso-potential lines

## EEG is a Measure of Neuronal Synchrony (1)


- potential produced by a single neuron is too small to be measured by an electrode on the scalp
- since the EEG-important neurons are vertically arranged, they can summate
- hence EEG is generated by a huge population of neurons:

cortical area required to produce a potential measured by a scalp electrode:  $\sim 6 \text{ cm}^2$ 

 $\sim 10^5$  neurons per 0.008 cm<sup>2</sup>

so the measured signal is from  $\sim 10^8$  neurons !!!

## EEG is a Measure of Neuronal Synchrony (2)




Time averaged potential P<sub>t</sub>


#### if *m* dipoles oscillate in synchrony:

$$P_t \approx m$$

if *m* dipoles oscillate non-synchronously:



## **EEG and Cortical Anatomy**




vertical-parallel arrangement of cortical pyramidal neurons is not always present (e.g. amygdala)

dipoles will tend to cancel each other out and hence the potential will approach zero for a large enough neuronal population

gyral surfaces contribute most !!

## What does the EEG measure?




current flow through wires vs. aqueous solutions:

wires

- negatively charged
- electrons move towards the positive charges

aqueous solutions - the ions do the flowing

## How do we know that Brain Generates the EEG?



- current flow follows the path of least resistance
- thus expect that current over skull defect will be greater than in a person with no skull defect

#### Adrian and Matthews (1934) "alpha-EEG"

Observation 1: alpha current density above skull defect is higher than in person with no defect (if alpha is generated outside the skull, this would not be observed)

Observation 2: less alpha current flows into frontal regions (if alpha is generated by eye muscles, this effect would not occur)


## Non-invasive measurements of EEG

## International 10/20 System:

- consistent naming
- reproducible placement



**Reference points:** 



	nasion (between forehead and nose)
	inion (bump at the back of the skull)
	left + right pre-auricular points (between lower jaw & cheek bone)
location:	F_rontal, T_emporal, P_arietal, O_ccipital, C_entral, z (central line)
numbers:	even: right side, odd: left side

## Non-invasive measurements of EEG

prepare scalp area (light abrasion to reduce impedance due to dead skin cells) place electrodes on scalp with a conductive gel or paste

place ground electrode (electrically zero?) and choose reference electrode(s) connect each electrode to one input of **differential amplifier** connect reference electrode(s) to other input of differential amplifier Amplification: x  $10^3$ - $10^5$ ; EEG has small amplitude: ~  $10-50 \mu$ V (AP: 60-100 mV) anti-aliasing filtering A-D conversion (typically 16 bit) low-pass: 35 - 70 Hz high-pass: 0.5 - 1 Hz sampling rate: 200 Hz (20 kHz) digital storage

**T**6

02

Common Refere

representation on screen

## Representation of EEG: the Choice of a Montage

EEG signal represents a difference between the voltages at two electrodes

#### Bipolar montage

difference between two adjacent electrodes (e.g. Fp1-F3)

#### **Referential montage**

difference between a certain electrode and a designated reference electrode

- midline positions (they do not amplify the signal in one hemisphere vs. the other)
- linked ears (physical or mathematical average of electrodes at both earlobes or mastoids)

#### Average reference montage

outputs of all of amplifiers are summed and averaged averaged signal is used as the common reference for each channel

#### Laplacian montage

difference between an electrode and a weighted average of the surrounding electrodes

#### awake, resting; normal posterior alpha rhythm disappears with eye opening (\*); high frequency activity after eye opening is muscle artifact; anterior-posterior bipolar montage; 0.5 - 70 Hz

- Man Marken Ma -----
- mound the property of the second second the second se F4-C4
- and marked and a second and a second and the second
- C3-P3
- when the second war war and the second war and the
- Mar mar and marked and the second an
- Marken Marke
- F8-T4
- 1 sec Fp2-F8
- month and the second se
- Marken M
- and many and a second and the second monor war have a second when the second when the second when the second se F7-T3

## An EEG Example

man man and the first of the second s

Market and a part of the stand of the stand

man man and the second s


## State-dependent Rhythms of the EEG


EEG power spectral density:

- frequency dependent amplitude behavior
- looks like 1/f noise
- long-term correlations

name f [Hz]	Α [μV]	Loc.
δ 0-4	20-200	variable
<b>9 4 - 8</b>	5-100	front., temp.
α 8-13	5-100	occ., par.
β 13 - 30	2-20	front.
γ 30 - 70	< 10	variable
(higher ??)		

functional integration: from short ( $\gamma$ ) to long-range interactions ( $\delta$ )






## Artifacts look more or less like EEG Signals

#### technical artifacts

- power line (incl. higher harmonics)
- broken cable, dirty connectors
- different electrode materials
- poor grounding

#### biological artifacts

- eye-induced (blink, movement)
- cardiac-induced (heart beat, pulse)
- muscle-induced (movement, chewing)
- sweating



## Artifacts look more or less like EEG Signals

CZ-C3	mm and here we have been and the second of t
СЗ-ТЗ	
[3-SP1	
SP1-SP2	wither where we
SP2-T4	Arren / "Monorman and the second and the
r4-C4	war war and and and the state of the state o
C4-CZ	man for an and the second seco
EOG1-EOG2	m for the second war and the sec
FP1-F7	min how we have a second with the second of
F7-T3	The back of the ba
[5-01	
FP2-F8	many white the second water and the second white the second with the second wi
-8-T4	have been been been been been been been be
F4-T6	manual of monorman and a service of the service of
F6-O2	and the second and th
FP1-F3	many manufacture and the second and
-3-C3	my for an and the second and the sec
C3-P3	un for the hard the part of th
P3-01	Manager and
FP2-F4	multimenter and the second of
-4-C4	way for the second when the se
C4-P4	and from the second and the second a
P4-02	
EKG	
	Eve Blink C3 pulse artifact C3
comment	Eye Blink C3 pulse artifact C3 pulse artifact C3 pulse artifact C3 pulse artifact Eyes left Eyes right 100 uV   1 sec

## Applications of EEG


well-established uses of EEG levels of consciousness (e.g. coma, delirium) epilepsy (e.g. diagnosis, epilepsy surgery) sleep (e.g. sleep disorders) other neurologic diseases, depth of anesthesia brain death, intensive care

areas of research

neuroscience, cognitive science, cognitive psychology, psychophysiological research brain-computer interfaces


many, many abuses

## EEG Examples: Sleep




Stage II sleep. K-complex (\*); Sleep spindles (\*\*). Anterior-posterior bipolar montage.

## EEG Examples: Primary Generalized Epilepsy




Burst of generalized 3 Hz spike and wave activity (\*). Anterior-posterior bipolar montage.

## EEG Examples: Primary Generalized Epilepsy




Transition to tonic-clonic seizure. Referential montage; reference = linked ears

## EEG Examples: Temporal Lobe Epilepsy



## EEG Examples: Temporal Lobe Epilepsy



## Summary

EEG is dominated by activity at the cortical surface (i.e., the gyrus) only about on-third of cortical surface

EEG cannot "see" deep into the brain spontaneous activity in, e.g., mesial temporal regions, interhemispheric frontal lobe structures, thalamus is NOT apparent on scalp EEG

Note that output of cortical neurons is in direction away from surface

suggests that is more sensitive to afferent (sensory) inputs to cortex (not proven)

Rhythmic nature of EEG unexplained !!!



 best spatial-temporal resolution (compared to all other measurements of cerebral activity)

temporal resolution: 1 msec (1000 Hz) or better spatial resolution: 10 microns (theoretically!)

- can measure while behavioral activities are on-going